

EN **RIALTO®**  
SURGICAL TECHNIQUE


A modified Kempe incision is recommended, as it ensures a better vascularisation of the myocutaneous flap and greater extension, thus helping the skin to seal up.

Draw the midline, from the hairline (widow's peak) to 1 cm above theinion (external occipital protuberance). The incision shall run parallel to the midline, 3 cm from it, medially. If the keyhole needs to be exposed, the incision may be cut 2 cm longer along the hairline, on the same side. Split the incision into two similarly sized flaps, closer to the ear. Spread out at the bottom, surrounding the pinna at the front, as far as just before the tragus.


Holes are usually made in the frontal-basal area (keyhole), at the base of the temple, in front of the tragus, and in the posterior portion, superomedially to theinion (HOLES IN RED). Additional holes along the upper edge (HOLES IN BLUE),


Make the widest possible craniotomy, avoiding the frontal sinus. The bone flap must include the frontal, parietal and temporal-occipital braincase and, as currently suggested by the Brain Trauma Foundation guidelines, the size must be at least 12x15 cm (surface >140 cm<sup>2</sup>).

**5 Cranial RIALTO brackets** are fitted in, starting from the skull and securing the part marked as C (Cranial). **Make a few guiding holes, approximately 3mm deep, where the 1mm tip needs to be inserted.**


Make sure the brackets are fitted in the right direction (with the letter C looking up). (FIGURE 4, Highlight the letter "C" better on the drawing).

The **brackets** are fitted in at the front, top and back of the edge of the craniotomy, not at the bottom (temporal bone). Once the 5 brackets have been fitted onto the skull, **gently and slowly mould them, without twisting**, using the pliers provided, until the bone flap can rest on the shelves of the brackets.


Once the brackets have been moulded so that the bone flap can steadily rest on the shelves, secure them with 2 screws each.

A or B) Different positions, depending on the thickness of the operculum

To help the soft tissues move closer to each other, usually you just need to place two strong, full-thickness crisscross pulling points at the corners of the front and back flaps.

Pln some selected cases, it may be helpful to detach the galea from the contra-lateral skull and/or make lengthwise and crosswise incisions in the over-turned flap of the galea and dilate it by hand.

Once decompression is no longer needed, the operculum **may be left in place**. If you decide to reposition it, reopen the flap of soft tissues, remove the flap of soft tissues, remove the **Cranial RIALTO® brackets** and secure the operculum, as in the traditional method.

REFERENCES

1) Carney N, Totten AM, O'Reilly C, Ullman JS, Hawryluk GW, Bell MJ, et al. *Guidelines for the Management of Severe Traumatic Brain Injury, Fourth Edition. Neurosurgery. 2017;80:6-15.*

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